



• Product Code

UTG-0415 Triple Motion Sieve Shaker with Time Adjustment, Motorized,
220-240 V 50-60 Hz

The UTG-0415 Triple Motion Sieve Shaker has a versatile design which allows the user to work with large batch sizes, large particle sizes or sieves with various sizes are need to be shaken.

This sieve shaker features a unique combination of jarring and orbital action which provides the most effective shaking action. Upper and lower cross heads are adjustable to be used for different sieve sizes. The shaker includes a 30 minute timer and a continuous-stop-timer operation switch.

• Technical Specifications

Sieve Capacity	10x of 200 mm (8") sieves plus pan and cover 8x of 250 mm (10") sieves plus pan and cover 6x of 300 mm (12") sieves plus pan and cover
Max. Sample Weight	Up to 4500 g depending on sieve size
Orbital Action	327 oscillations per minute (approx.)
Jarring Action	40 vertical blows per minute
Dimensions	540x370x1015 mm
Weight (approx.)	75 kg
Power	250 W

