



BPS INSTRUMENT

TENSILE TESTING MACHINE

• Product Code

UTG-0418 Motorized Sieve Shaker, 3", 5" and 8" Sieves, 220-240 V 50-60 Hz
UTG-0418/110 Motorized Sieve Shaker, 3", 5" and 8" Sieves, 110 V 60 Hz

• Standards

The UTG-0418 Motorized Sieve Shaker imparts a uniform shaking action to a range of sieve dimensions. The shaker can handle up to 10 X 8", 12 X 5" and 16 X 3" full-height sieves and, 18 X 8" half height sieves. Mounting is required.

Timer
30 minute, analogue
Dimensions
380x380x1145 mm
Weight (approx.)
30 kg
Power
1/4 HP

